

everest biolabs

Automate the tedious.
Focus on discovery.

Apex[®]

Easily isolate pure, functional
EVs with high yield



Also available
in plate and
mini formats




Plate and column
formats



SEC tuned for
purity and yield



Reproducible
< 10% CV



0.1 - 1.0 mL
sample volumes



Apex for sharper peaks. Cleaner insights.


Apex 4B and 6B columns & plates

High performance isolation

- Reproducible performance < 10% CV
- Sample volume 0.5 - 1.0 mL
- Available with SEC and multi-mode resins
- Validated in culture, plasma, urine, serum & CSF
- Compatible with Ascent and Summit

Column Type	Purity	Yield	Pore-Size
Apex 6B	1.0	1.0	20 nm
Apex 4B	5.0	0.8	35 nm
Competitor	5.0	0.3	70 nm


EV profiles from 0.5 mL plasma, as measured by Atlas ELISA, after SEC with agarose gels of 70 nm, 35 nm, and 20 nm pore sizes. Definitions: Yield = EVs in fractions 1-3, normalized to Apex 6B Purity = Yield/HSA, normalized to Apex 6B



Apex MM, more than just SEC

50% more EVs and 2x purity


SEC combined with multi-mode (MM) resin significantly enhance EV purity by removing more lipoproteins than SEC alone.



APEX SPECIFICATIONS			
Column format	Column	Plate	Mini
Column type	MM / 6B / 4B	6B / 4B	MM
Input sample volume	0.5 - 1.0 mL	0.5 - 1.0 mL	0.1 - 0.2 mL
Column reproducibility	10% CV		
Sample types	Plasma, serum, urine, CSF, cell culture media		
Resin types	4 % or 6 % cross linked agarose beads, cross-linked agarose plus multi-mode		
Column volume	8.75-9.0 mL	8.75-9.0 mL	3.5 mL

Apex 4B

Optimized for high purity and lowest protein co-elution




Apex 6B

Optimized for highest EV yield



Apex MM

Optimized for highest purity



Apex Mini

Optimized for small sample volumes (0.1 - 0.2 mL)



See Apex performance data (NTA/TEM/ONI/...)

